

# SKB1A02

## FOLDABLE STRETCHER

Registration code: SKA-1B



# SKB1A02

## FEATURES AND HIGHLIGHTS

### Technical parameters

▪ Product size (LxWxH)	2200x530x140mm
▪ Folded size (LxWxH)	1100x180x130mm
▪ Loading bearing	159kg
▪ Net weight	6kg

#### Design

Scientific design, ergonomics, simple structure, simple operation, reliable and secure. Widely used in the home for the old people, school security, hospitals, hotels etc for emergency evacuation.




#### Leather

Made of high-strength oxford leather, fireproof, yellowing resistant, mildew resistant, tortuous, wear resistant 100,000 times.



#### Material

High strength aluminum alloy, foldable stretcher, compact and lightweight, easy to carry.



#### Function

Foldable stretcher, 2-fold. easy for carrying and storage.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue