



Grace Care series

D6w

ELECTRIC BED


Registration code: SKD-C



SAIKANG / Medical

Healthier
your life

D6w Electric bed



Technical parameters

▪ Size (LxW):	2245x1000 mm
▪ Bed height adjustment:	435-745 mm
▪ Mattress platform:	1925x830 mm
▪ Safe working load:	250 kg
▪ Back-rest adjustment:	0-75° (±10°)
▪ Knee-rest adjustment:	0-35° (±10°)

Technical configuration

<input checked="" type="checkbox"/> Linear motor	3 pcs.
<input checked="" type="checkbox"/> Hand controller	1 set
<input checked="" type="checkbox"/> 5" hardshell castors with break	4 pcs.
<input checked="" type="checkbox"/> Foldable side rail with ProErgo-Handle™	1 set
<input checked="" type="checkbox"/> ABS bed ends with safe lock	1 set
<input checked="" type="checkbox"/> IV pole	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs.
<input checked="" type="checkbox"/> Drainage hook	2 pcs.

Optional

<input type="checkbox"/> Manual CPR
<input type="checkbox"/> Back-up battery
<input type="checkbox"/> X-ray translucent mattress platform
<input type="checkbox"/> Directional wheels function
<input type="checkbox"/> Monkey pole
<input type="checkbox"/> Aluminium balkan frame



SAIKANG Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 +86 18021231901 · 24/7
 sales@chinasaikang.com
 www.saikangmedical.com



Company video



General catalogue