



A100-4

Electro-Hydraulic
Operating Table

SKL-A



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

DIMENSION PARAMETER

External Length	2100mm
External Width	500mm
Table Height Min.	800mm
Table Height Max.	1000mm

TECHNICAL PARAMETERS

Reverse Trendelenburg/Trendelenburg	22°/25°
Lateral Tilt Right/Left	20°/20°
Head Rest (Up/Down)	45°/90°
Foot Plate (Up/Down)	15°/90°
Foot Plate Horizontal Positioning	90°
Back Plate (Up/Down)	75°/20°
Kidney Bridge Elevation	110mm
Horizontal Sliding	300mm

STANDARD ACCESSORIES

Bed Frame	1pc
Mattress	1set
Anesthesia Screen Frame	1pc
Shoulder Holder with Cushion	2pcs
Arm Holder with Cushion	2pcs
Leg Support with Cushion	2pcs
Foot Support with Cushion	2pcs
Kidney Bridge Handle	1pc

QUALITY CONTROL



PROFESSIONAL FUNCTIONS

Lowest Position of
Table Top

750mm



Highest Position of
Table Top

1050mm



(Optional) Head Rest



Reversed Trendelenburg

0°-22°



Trendelenburg

0°-25°



(Optional) Lateral Arm
Support



Right Lateral Tilt

0°-20°



Left Lateral Tilt

0°-20°



(Optional) Waste Basin
For Gynecology



Anorectum
Operation Position



E.N.T Operation Position



(Optional) Orthopedics
Traction Frame



RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



Armrest



Anesthesia Frame



I.V Pole



Leg Crutches



Body Supporter





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.