



# GB8e

## ORTHOPEDIC TRACTION ELECTRIC BED

Registration code: SKD-C




**SAIKANG** Medical

Healthier  
your life

# GB8e

## Orthopedic traction electric bed



### Technical parameters

▪ Size (LxW):	2250x1060 mm
▪ Height adjustment of use base:	435-725 mm
▪ Mattress platform:	1935x900 mm
▪ Safe working load:	250 kg
▪ Back-rest adjustment:	0-75° (±10°)
▪ Knee-rest adjustment:	0-35° (±10°)
▪ Reverse trendelenburg:	0-12°
▪ Trendelenburg:	0-12°

### Technical configuration

<input checked="" type="checkbox"/> Linear motor	3 pcs.
<input checked="" type="checkbox"/> Hand controller	1 set
<input checked="" type="checkbox"/> Angle indicator	4 pcs.
<input checked="" type="checkbox"/> Double side castors (5") with central locking system	4 pcs.
<input checked="" type="checkbox"/> Central locking pedal	2 pcs.
<input checked="" type="checkbox"/> ABS side rails with extra height	1 set
<input checked="" type="checkbox"/> ABS bed ends with safe lock	1 set
<input checked="" type="checkbox"/> Aluminium balkan frame	1 set
<input checked="" type="checkbox"/> IV pole	1 set
<input checked="" type="checkbox"/> IV pole prevision	4 pcs.
<input checked="" type="checkbox"/> Drainage hook	2 pcs.



#### Virgin ABS

Environmental plastic can be sterilized in 100°C, tensile upto 30 MPa.



#### Bed surface

Stamping forming platform, smooth surface no welding scar, bearing 250 kg load. Option ABS platform and resin X-ray translucent platform.



#### Steel and painting

Baosteel from Fortune Global 500 11 Process epoxy painting, ASTM testing anti-bacterial, paint thickness 0.12 mm, brightness 60°, paint can resist 50 kg impact.



**SAIKANG** Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China  
 📞 +86 18021231901 · 🕒 24/7  
 ✉ sales@chinasaikang.com  
 🌐 www.saikangmedical.com



Company video



General catalogue