



SKR-IV588

INFUSION TROLLEY

Registration code: SKM-B

SKR-IV588

FEATURES AND HIGHLIGHTS

Technical parameters

- External size (LxWxH) 588x480x891 mm

Technical configuration

- Cabinet 1pc
- Drawers 2pcs
- Casters 4pcs
- Trash can 2pcs


Material

Made of the E1 Class HPL board and aluminum alloy frame, resistance of water, acid, alkali, fire and easy clean, not easily polluted by reagents, blood, drugs etc., satisfy the daily disinfectant and infection control with long life.



Guardrail

Aluminum alloy guardrail on three side, anti-corrosion and anti-rust



Drawer

Two drawers and 1 cabinet made of high quality HPL board, with high strength aluminum handles



I.V Pole

Height adjustable I.V pole



Infusion Rack

Fixed double-row 5-hook infusion rack



Caster

TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts
Pass Dynamic test: Bearing 120kg run 30KM and pass obstacles 500time



Trash Can

2 ABS trash can without lid





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai,
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China
 📞 +86 18021231901 · 🕒 24/7
 ✉️ sales@chinasaikang.com
 🌐 www.saikangmedical.com