



SKR054-ET

Emergency Trolley



SAIKANG / Medical



LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098

QUALITY CONTROL

TABLE TOP

With armrests on both sides, professional sharp container, can be placed left and right, surface with transparent soft glass



CHASSIS

Smaller than the table top size, saving the bottom space, reduce the collision during the movement, more convenient

DRAWER

Drawer insize big & deeper!

- *Two small drawer panel height: 80mm, Inner size: 430*335*68mm
- *Two middle drawer panel height: 120mm, Inner size: 430*335*110mm
- *One big drawer panel height: 240mm, Inner size: 430*335*220mm



CENTRALIZED LOCK

Inner with partitions can be well organized

POWER SOCKET

The power socket is movable and detachable, easy to use, located above the CPR board to convenient the use of oxygen cylinders



OXYGEN BOTTLE HOLDER

Inserted retractable oxygen bottle holder, belts fixed on trolley main body, space saving

DEFIBRILLATOR

The plate surface is made of virgin ABS material. The position of the I,V pole and the defibrillator can be exchanged for easy use.




TRASH CAN

Two dusk baskets are divided into two colors, convenient for medical staff to classify different types of garbage.

CASTERS

4 inch luxurious Noiseless casters, and two with brake
High strength polyurethane material, no worn out after running 30KM, save anti-winding hard shell, united forming without bolts, portable and flexible



DRAWER CARD HOLDER

Card insert from downside, easy keep clean.





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.