




SKB037(C)


Patient Trolley



SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE

Stock Code: 870098



Back-rest Knee-rest Leg Plate Flating Height Adjustment

MAIN TECHNICAL PARAMETERS

| | |
|-------------------|----------------------|
| Size | 1960*700*(620-880)mm |
| Back-rest | 0-75° |
| Knee-rest | 0-35° |
| Leg Plate Flating | 0-5° |

TECHNICAL CONFIGURATION

| | |
|------------------------------------|------|
| Aluminum alloy side rail Side Rail | 1pc |
| 5" Covered Castors | 4pcs |
| IV Pole | 1set |
| Handle Crank | 3pcs |
| Debris Plate | 1pc |
| Oxygen Bottle Holder | 1pc |

QUALITY CONTROL

BED FRAME

The whole frame used high quality S.S. with the mirror surface polishing



BED SURFACE

The stretcher surface made by phenolic compact board, much thinner and easier to X-ray



HANDLE CRANK

The bed with three function, controlled by manual crank



BUMPER STRIP

The bumper strip add extra 6pcs rivets on each side to make sure the fastness.



GUARDRAIL

Aluminum alloy side rail, 5-staff, knocked-down type. Tight by screw, easy to change if it is broken.



RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



5" silence caster, advanced protection universal brake wheel, easy to move, simple to operate, reliable performance



Adjustable stainless steel I.V pole, outside diameter Φ 19*1, Inside diameter Φ 16*1, Load capacity: 5KG, Fixed device was added and ensure the I.V pole will be very stable with no shaking.



1pc Debris Plate



1pc Oxygen Bottle Holder





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.