



SKB041-3

Patient
Transportation
Trolley




SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098



Back-rest Knee-rest Height Adjustment



Reverse Trendelenburg Trendelenburg

MAIN TECHNICAL PARAMETERS

Size	2100×850×(630-930)mm
Back-rest	0-75°(±10°)
Knee-rest	0-45°
Reverse Trendelenburg	0-15°
Trendelenburg	0-15°
Safe Working Load	250kg

TECHNICAL CONFIGURATION

USA Hydraulic Pump	2pcs
Crank System	1 set
Back Rest Air Spring	2pcs
8" German Castors,Central Locking	4pcs
Emergency Push Handle	1 set
Directional Wheel	1 pc
Folding Side Rail	1 set
Φ 74mm Wall Bumper Wheel	1 pcs
IV Pole Previsions	2pcs
Mattress	1 pc
IV Pole	1 set
Central Locking Break Pedal	2sets
Directional Wheel Release Pedal	2pcs
X-Ray Translucent Bed Platform	1 set

QUALITY CONTROL

METAL MATERIAL

Baosteel from Fortune Global-500
Panasonic Robot ensure 360° full
smooth welding

500 | 360°

PLASTIC

Environmental plastic can
be sterilized in 100°C,
tensile upto 30MPa

100. | 30

STEELS | PAINTING

1 | Process epoxy painting,
ASTM testing anti-bacterial, paint
thickness 0.12mm, brightness
60°, paint can resist 50kg impact

0.12... | 60° | 50.

BRAKE


Break spring keep qualified
after 10,000times
uninteruped break

HYDRAULIC PUMPS

Height adjustment,
trendelenburg and
revised trendelenburg
through two USA
hydraulic pumps

CENTER WHEEL

Can be engaged to facilitate
direction control, enabling
users to navigate the
stretcher easily in hectic and
dynamic environments.



RICH PARTS

LEADING BRAND OF
CHINA MEDICAL FURNITURE



Hydraulic treadle for rise



- ① Leg section lows down to make head highter (Than Leg Section)
- ② Back section lows down to make leg section highter (Than head Section)
- ③ Step on the middle of the pedal, make the bed base fall down steadily



Four double side casters of 6 inch diameter; simultaneously lockable by the step pedal at both sides.
TPR tire no worn out after running 30KM, save anti-winding hard shell, united forming without bolts



Waterproof PU maffress, convenient for transfer patient



Foldable aluminium alloy guardrail



Bumper wheel detachable, max. protect the bed and patient during transportation



Full SS IV pole, bearing 15kg





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03-England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.