



SKE180

Electric wheelchair




SAIKANG / Medical

LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

QUALITY CONTROL

DESIGN

Scientific design, ergonomics, providing a comfortable experience.
Foldable type for easy transport



LITHIUM BATTERY

High-performance lithium battery, detachable, safe and reliable, range up to 20 kilometers,



REMOVABLE CONTROLLER PANEL

With five-speed adjustment, horn button, and LED light display power, Joystick control direction, easy to operate



FOOT PANEL

High quality plastic foot panel, removeable and foldable



WHEEL

Solid rubber universal wheels, high quality and reliable.
Different sizes wheels are available, meet the request of different patients.



HANDLE BRAKE

Easy to operate, safe and reliable performance



CLUTCH

Each motor with a clutch, The wheelchair can be switched between electric and manual by operating the clutch





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-1802123190 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

Uk Healthward International Limited
Owner(s) Address DEPT: 114,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Contact



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.