



SKH041(1)

Stainless Steel
I.V. pole



SAIKANG / Medical


LEADING BRAND OF CHINA MEDICAL FURNITURE
Stock Code: 870098

TECHNICAL PARAMETERS

Size: $\Phi 25 \times 1350-2250$ mm

TECHNICAL CONFIGURATION

2" Pluger Wheel	5pcs
Plastic Pentagram Caster	1pc



QUALITY CONTROL

DESIGN

Movable stainless steel I.V. pole with four hooks



HEIGHT ADJUSTMENT

Height adjustable with side knob



INFUSION HOOK

Four collapsible infusion hooks, single load-bearing 1 kg



WHEELS

5-star good quality ABS base, with castors for flexible and stable movement





INNOVATION OF SAIKANG

LEADING BRAND OF CHINA MEDICAL FURNITURE



Design Innovation



Research Innovation



Manufacturing Innovation



Maintenance Innovation



Service Innovation

JIANGSU SAIKANG MEDICAL EQUIPMENT CO.,LTD.

24hrs: 0086-18021231901 
www.chinasaikang.com

SAIKANG (China Headquarter)

No.35 Lehong Road,Zhangjiagang City,
Jiangsu Province,China.
TEL: +86-512-58525850
FAX: +86-512-58520713
E-mail: inquiry@jssaikang.com

SAIKANG (USA Branch)

Saikang Medical Technology, Inc
17800 Castleton St Ste 180,
City of Industry, Ca 91748
Tel: +1 949648-2776

SAIKANG (UK Branch)

UK Healthward International Limited
Owner(s) Address-DEPT 111,
196 High Road Wood Green,
London, N22 8HH
United Kingdom
TEL: +44 2032399738

SAIKANG (East Africa Branch)

Saikang Medical Solutions (East Africa) Ltd.
Viraj Complex Unit 3 Mombasa Road,
PO Box 14262-00800, Nairobi.
TEL: +254 733-704 690

SAIKANG (Mid East Branch)

Saikang Medical Supplies (Mid East) LLC
X22-S03 England Cluster,
International City, Dubai UAE
TEL: +20 128 8962497



Catalogue



Website

As the circumstance of printing, the photos might have minor difference on color, tactile from real products.
As the update of design and improvement, parameters and outlook of products might change, thanks for understanding.