


SKH046-2

OVERBED TABLE



SKH046-2

FEATURES AND HIGHLIGHTS



Technical parameters

- External size (LxWxH) 890×440×(765-1050)mm

Technical configuration

- Air Spring for Height Adjustment 1pc
- ABS Table Top with Cup Socket 1pc


Application

Movable dining table, nursing bedside table, can be used with various hospital beds.



Table Top

ABS material table top, easy to clean, endurable to chemical cleaning. With cup holder, and edge rising up for protecting the falling. Round corner design, anti-collision and durable.



Height Adjustment

The gas spring system help the table cover moving up and down.



Bottom Frame

The bottom frame is sprayed with carbon steel and the center of gravity is stable.



Castors

4 castors made by rubber and could revolve around, easy to install, more stable and firmer..





CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue