



# SKR-CT232

## CLINICAL TROLLEY

Registration code: SKM-B



**SAIKANG** / Medical

Healthier  
your life

# SKR-CT232

## FEATURES AND HIGHLIGHTS

### Technical parameters

- Size (LxWxH) 780×475×1000 mm

### Technical configuration

<input checked="" type="checkbox"/> File Box	1pc
<input checked="" type="checkbox"/> Plastic Basket	2pcs
<input checked="" type="checkbox"/> Sharps Boxes	2pcs
<input checked="" type="checkbox"/> Trash Can	2pcs
<input checked="" type="checkbox"/> Castor	4pcs



#### Table Top

Guardrails on top of trolley is made of stainless steel, table top with soft glass, avoid possible movement and fall down during movement of trolley.



#### Material

Steel and aluminum structure inside, ABS materials outside, easy to clean.



#### Push Handle

With push handle on both side, convenient for pushing easily.




#### Drawer

1 layer small drawer, 2-3 layer middle drawer, can put different sizes of medicine boxes and other medical materials, and equipped with central control lock.



#### Storage Shelf

Storage box which is designed fixed on the frame, large volume, easy to store and take out items, easy to clean.



#### Drawer Inside Dividers

With abs inserts inside, which can be separated freely as required.



#### Drawer Card Holder

Card insert from downside, easy keep clean.



#### File Box

Side attached file box, for storage of paper work.



#### Plastic Basket

With two plastic basket, there are two 2-liter Sharps boxes in the basket.



#### Trash Can

Two trash cans in different color with lid, controlled by foot switch.



#### Castor

4 nylon wheels Ø100 mm with non-marking tires and plain bearings, two with brakes.





## CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



## WORLDWIDE BRANCHES

### USA

**Saikang Medical Technology, Inc.**  
17800 Castleton St Ste 180,  
City of Industry, 91748, California.  
Tel.: +1 949648 2776

### United Kingdom

**UK Healthward International Ltd.**  
Dept. 111, 196 High Road Wood  
Green, N22 8HH, London.  
Tel.: +44 2032399738

### Russia

**Saikang Medical Russia LLC**  
Truda 174,  
454080 Chelyabinsk.  
Tel.: +7 (351) 214 5566 / 8811 / 1516

### East Africa

**Saikang Medical Solutions Ltd.**  
Viraj Complex Unit 3, Mombasa  
Road, P.O. Box 14262-00800, Nairobi.  
Kenya.  
Tel.: +254 733 704 690

### Middle East

**Saikang Medical Supplies LLC**  
X22-S03 England Cluster,  
International City, Dubai.  
United Arab Emirates.  
Tel.: +20 128 8962497

# SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.  
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎️ +86 18021231901 · 🕒 24/7

✉️ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue