



A105

OBSTETRIC TABLE

Registration code: SKF-C




SAIKANG / Medical

Healthier
your life

A105

FEATURES AND HIGHLIGHTS



Technical parameters

External Size (LxWxH)	1900x600x700-900 mm
Back Plate Up/Down	50°/3°
Trendelenburg/Reverse Trendelenburg	10°/20°
Leg Board Outwards	90°
Working Area	600x600 mm
Patient Weight Capacity	185kg


Technical configuration

<input checked="" type="checkbox"/> Bed Frame	1pc
<input checked="" type="checkbox"/> Mattress	1pc
<input checked="" type="checkbox"/> Hand Bars	2pcs
<input checked="" type="checkbox"/> 2 inch casters with brake	2pcs
<input checked="" type="checkbox"/> Arm Support with Cushion	2pcs
<input checked="" type="checkbox"/> Leg Support with Cushion	2pcs
<input checked="" type="checkbox"/> Anesthesia Screen Frame	1pc
<input checked="" type="checkbox"/> Stainless Steel Basin	1pc
<input checked="" type="checkbox"/> Handwheel	4pcs
<input checked="" type="checkbox"/> Oil Pot (With Oil)	1pc



Material


The base seat and upright column covering are made of high quality 304 stainless steel, which is anti-corrosive and can be easily cleaned.
Optional: Powder Coated Steel Base.

**Function**

The main bed and auxiliary bed are connected with each other, auxiliary table is easy for rotating and dismantling.

**Rotary Brake Device**

Rotary brake device which is easy for moving obstetric table.

**Foot Step Hydraulic Elevating**

Foot step hydraulic elevating, Imported Y type sealing ring and is durable.

**Mattress**

The mattress adopts high density memory sponge and is once formed, with no seam.
Optional: Detachable mattress.

**Stainless Steel Sewage Basin**

Equipped with stainless steel sewage basin and other various accessories.

**Leg Support with Cushion**

Double using leg support, can be used as boatstretcher.

**Anesthesia Screen Frame****Arm Support with Cushion****Hand Bars****Handwheel**



CONSTANT INNOVATION



Design



Research



Manufacturing



Maintenance



Service



WORLDWIDE BRANCHES

USA

Saikang Medical Technology, Inc.
17800 Castleton St Ste 180,
City of Industry, 91748, California.
Tel.: +1 949648 2776

United Kingdom

UK Healthward International Ltd.
Dept. 111, 196 High Road Wood
Green, N22 8HH, London.
Tel.: +44 2032399738

Russia

Saikang Medical Russia LLC
Truda 174,
454080 Chelyabinsk.
Tel.: +7 (351) 214 5566 / 8811 / 1516

East Africa

Saikang Medical Solutions Ltd.
Viraj Complex Unit 3, Mombasa
Road, P.O. Box 14262-00800, Nairobi.
Kenya.
Tel.: +254 733 704 690

Middle East

Saikang Medical Supplies LLC
X22-S03 England Cluster,
International City, Dubai.
United Arab Emirates.
Tel.: +20 128 8962497

SAIKANG / Medical

Jiangsu Saikang Medical Equipment Co., Ltd.
STOCK CODE: 870098

📍 No. 35, Lehong Road · Zhangjiagang (Jiangsu) · China

☎ +86 18021231901 · 📞 24/7

✉ sales@chinasaikang.com

🌐 www.saikangmedical.com



Company video



General catalogue